

Keep all your audiology equipment  
clean and patient-ready

# MEDTROL'S AUDIOTECH® CLEANSING SURFACE WIPES





**PROTECT  
YOUR EQUIPMENT,  
INSTRUMENTATION AND  
YOUR PROFESSIONAL IMAGE  
WITH AUDIOTECH**

## Product Description:

AUDIOTECH contains a non-alcohol, pre-saturated (quaternary ammonium chloride) wipe designed to clean and deodorize all surfaces and is specifically recommended for audiology equipment and instrumentation. The no-rinse, quick drying solution reduces oil build-up and is compatible with plastics, rubber, metals, clothing and all fabrics.

**Suggested Areas of Use:** This fresh scented sanitary cleaner and deodorizer can be used on plastics, rubber, hearing aids, testing devices, assisted listening devices, ear molds, specula, otoscopes, cell phones, ear phones and all non-invasive audiology instrumentation.

Product	Name	Catalog #	Item	Unit Size	Case Size
	AUDIOTECH	AT217	6"x 7" wipe	160 wipes per canister	12 canisters / case (1,920 wipes)
	AUDIOTECH Mini	AT219	4"x 6.5" wipe	36 wipes per canister	24 canisters / case (864 wipes)

*Private labeling available  
for this product.*



HEALTHCARE DIVISION

7157 N. Austin Ave. • Niles, IL 60714  
Ph: 800-647-7180 • Fx: 800-255-3027  
sales@medtrol.com • www.medtrol.com

© 2008 Medtrol, Inc.

**CONTACT MEDTROL TODAY AT 800-647-7180 TO LEARN MORE!**